


## LabTIE DISPENSERS

Simple. Fast. Reliable.

### Product Info



## Introduction

LabTIE Dispensers are used to quickly dispense Seeds, Beads, Microbeads and/or Powders in industry standard plates, tubes, vials or various trays. Due to the modular design of exchangeable Mesh Plates the Dispenser can be used to dispense a variety of sizes, amounts and volumes.

- Seeds 0.2mm – 9.0mm
- Beads 0.3mm – 6.0mm
- Powders 1uL/1mg – 1500uL/1.5gr

Each Dispenser is designed, produced and tested by our own Tech Division and custom made to the requirements of our customer. Our Dispensers are used all over the world in research laboratories, universities, pharmaceutical– and biochemical laboratories, institutes and by seed breeding companies.

## Advantages:


- Dispensing in only 30 seconds
- Accurate dispensing, <1% error rate
- Easy cleaning
- Exchangeable volume or seed Meshplates plates, “Different volume or Different seed variety? Just slide in a different volume plate”
- Custom made volume and seed plates available
- Dispenser can be customized to fit every ‘container’ where you want to dispense your solid object in
- Sterile dispensing



## Dimensions

All standard Dispenser Dimensions can be found in the PDF on [www.labtie.com](http://www.labtie.com) as the corresponding Dispenser product. When a Dispenser is Custom made, dimensions may vary, ask your LabTIE contact about the specifications.

## Product Composition



## Standard Plates, Tubes and Trays supported by LabTIE Dispensers

**Is your 'container' not listed? Dispenser can be customized to fit your 'container'**

### Microtiter and Cell Culture Plates



96-Well Microtiter, Cell Culture, Deep Well Plates and Plates with strips



384-Well Microtiter and Deep Well Plates



24-Well Greiner/Corning/ThermoF Plates    48-Well Greiner/Corning/ThermoF Plates

### Tubes



Snap cap tube 0.2mL-1.5mL



Screw cap tube 1.5-2.0mL  
Cryo Tube



1.5-2.0mL Autosample (GC)  
Vials

## Petri Dishes Round and Square



9cm and 10cm  
Agar/Blot paper




10cm and 12cm  
Agar/Blot paper

## Trays



Germination Tray  
Blot paper/Soil



108-Well Tray  
Soil



Tray Soil 35 x 22 cm

Any questions? Feel free to contact us:



**AXT PTY LTD**  
Authorised Distributor  
LabTIE  
Australia & New Zealand

1/3 Vuko Place  
Warriewood  
NSW 2102 Australia

+61 (0)2 9450 1359  
axt.com.au  
info@axt.com.au