

## Parallel synthesis the easy way

You need a flexible reaction block for manifold applications? You should try a VariSynth. With one instrument you can run different reactions in parallel, from 2 to 1000 ml.

The VariSynth System is available with 2 or 4 zones. Each can run a 1000 ml reaction or multiple smaller reactions up to 20x 2 ml with our special inlets.

Each zone of the VariSynth works independently. You can run up to four similar or different reactions in parallel without time-consuming preparation. Everything is operated by a large colored touch display which is very easy to use.

## Specifications

- -70 °C to +250 °C
- 2 or 4 zones
- 80 to 250 volt AC
- colored touch display
- 2 to 1000 ml volume
- 100 to 1500 rpm
- accessories    glassware  
                         stirrers  
                         sensors

## Contact

LAIX Technologies UG  
Am Fasanenhang 5  
52379 Langerwehe

Tel.: +49 (0) 2409 48798 07  
Fax: +49 (0) 2409 48798 08

Email: [info@laix-tech.de](mailto:info@laix-tech.de)  
Web: [www.laix-tech.de](http://www.laix-tech.de)



**Touch display**  
**Flexible RV size**  
**Suitable accessories**